

Assistance for those suffering from conditions leading to loss of dexterity or muscle strength. A call button YOU DO NOT HAVE TO PRESS.



Nurse Call

Versions available for operating established Nurse Call or for use with our Safelink radio pager

Assistance also for the carer by giving to two levels of urgency to the alarm

- Simply resting any part of your body against the front for more than 1 second will trigger an "ASSIST ME" alarm on the pager.
- Keep activating the area for more than 5 seconds and the alarm is elevated to "URGENT".



Or for more severe situations ask about using our sound triggered Companion monitor with the unique background sound rejection filter.

This is just one of the many distress sensing products available, which all communicate the alarm condition to the P137 pager to provide carers with extended vigilance.