



The P117 Autodialler is an add-on peripheral to most Alert-it Alarm Systems and is compatible with wired and radio systems.

When an alarm condition is detected, the dialler is triggered either directly in wired systems or via the radio (P117C). The dialler will then phone a series of programmed telephone numbers in sequence, delivering a short recorded audio message as to the situation. The dialling stops when the call is acknowledged by a carer

The dialler is able to warn of monitor problems (eg the need for a new battery) as well as user alarms

Most of the Alert-it Systems are set to give the user an audible warning for 20 seconds before the dialler is activated, so that false alarms can be intercepted

Alert-it Care Alarms are social aids designed and manufactured in accordance with 93/42/EEC as Class 1 Medical Devices. They are intended to improve, but not replace, the vigilance of carers to distressing side-effects of various diseases, such as Epilepsy and Dementia. They do not monitor vital physiological processes and should not be expected to diagnose any disease or predict the onset of any symptoms.