



ENURESIS SUPPORT



A unique moisture detector system with sensors embedded into thin cotton sheets. Detects moisture from urination



- Cotton sheet for user comfort and discretion
- Also reduces sweating, bed sores and false alarms
- Placed on a mattress or pillow
- Contains anti-allergenic silver coated wires
- Simple and easy to install
- Supplied with a detector unit and advanced pager
- Suitable for carers at home or professional care organisations.

Versions for direct
Nurse Call
connection available



Made in the UK



Stress free caring ...



MOISTURE DETECTOR

The Comf-it moisture detector system is ideal for detecting when a vulnerable person has either urinated or vomited. The system alerts the carer to a build up of moisture via a small detector unit. The alarm can either be transmitted through a failsafe radio link to an advanced pager with a range of 450m, or directly to a wired Nurse Call system.

The cotton sheets can be washed at normal temperatures once connecting cables are removed. The sheet size is 90 x 60cm with a 30cm side extension to tuck over the mattress.

RADIO DETECTOR UNIT



The small detector unit is battery powered with a rechargeable option. An integrated call button also allows the user to request assistance.

In the radio system the pager distinguishes between the different severities of each alarm and is powered by rechargeable batteries. For professional use a 32 channel version is available.

In the wired system different connection leads can be supplied to match any Nurse Call system.

WIRED DETECTOR UNIT



ADVANCED RADIO PAGER

For further information or to order the right system, please call us on 0845 217 9952, or visit our website at www.alert-it.co.uk