

Are you caring for someone with Dementia?

- Are you concerned to know if they have slipped outside when you were not watching?
- Are they at risk if they went into the kitchen or bathroom unescorted?
- Are they likely to wander around at night?

Over
10% less
than buying
the parts

The door switch can be disarmed with the secure key supplied



Two door activity switches (for Front and Back Doors)

- One motion sensor that can be used to:
- Detect entry into a high risk area (kitchen or bathroom)
 - Detect someone getting out of bed
 - Detect access to corridors or stair wells.

The Motion Sensor is battery powered and easily mounted in a discrete tamper-proof location. It is armed/disarmed by a small keyfob



SafeLink



The Safelink Pager system is supplied to Local Authorities, Hospitals, Residential Homes and Schools.

A small simple pager with 450m range to alert the carer to any risk situation, using high integrity radio

There are many other monitors compatible with the pager that can be added to the system if required. The P135 Alarm Badge, for instance, can detect falls and if someone wanders outside a safe area (eg the garden). All systems are covered by our 1 month risk-free trial period. We can also support those living alone by telephone based systems